

Metabolism

Measurement of bacterial metabolism

Risk assessment is valid for

Determination of metabolism

Analysis

The risk associated with this procedure will be minimized by following of this protocol.

List of Chemicals and Reagents

| Chemical/reagent | Link | Stock/ User Concentration | Special Waste Disposal?* |
|-------------------------|----------------------|---------------------------|--------------------------|
| Sulfuric acid | MSDS | 18M / 0.00125M | Y/N |
| Isotope ^{14}C | MSDS | 1mCi / 0.1 μCi | Y |

Secure Job Analysis

| Tasks | Unwanted Incidents | Precautions | Actions |
|--|-----------------------------------|---|---|
| Dilution of Isotope ^{14}C 1mCi | Inhalation and Skin contamination | Use double nitrile gloves and mouth masks. The work should be done in fume hood. | Wash the contaminated local with Plenty of water. Rinse eyes with water and Contact emergency if necessary. |
| Dilution of strong sulfuric acid 97% 18M | Burns on skin/eyes | Use nitrile gloves and Dilution of concentrated acid should be done in fume hood | Cool skin with cold water. Rinse eyes with water using eyewash spray head. Contact emergency if necessary. |
| Usage of 0.00125M sulfuric acid as mobile phase in HPLC system | Burns on skin/eyes | Use nitrile gloves | Cool skin with cold water. Rinse eyes with water using eyewash spray head. Contact emergency if necessary. |

*Check [SOP](#) (Standard operating procedure) for further details

HPLC

SHIMADZU High Performance Liquid Chromatography (HPLC)

Risk assessment is valid for

Operation of HPLC

Analysis

If this protocol is followed, there is **minimal risk** associated with use of this procedure.

List of Chemicals and Reagents

| Chemical/reagent | Link | Stock/ User Concentration | Special Waste Disposal |
|------------------|------|---------------------------|------------------------|
| None | | None | None |

Secure Job Analysis

| Tasks | Unwanted Incidents | Precautions | Actions |
|--------------------|--|--------------------------------|---|
| Loading of Samples | Spill of the isotope (^{14}C) mix from injector | Cover with foil during reading | Use wipe paper to wipe off spill use 70% ethanol as decontaminant and disinfectant. Sort the contaminated paper as radioactive trash. Cut the contaminated area of bench paper and Sort it as radioactive trash |

Centrifuge

Centrifuges

Risk assessment is valid for

Operation of centrifuges

Analysis

If this protocol is followed, there is **minimal risk** associated with use of this procedure.

List of Chemicals and Reagents

| Chemical/reagent | Link | Stock/ User Concentration | Special Waste Disposal |
|------------------|----------------------|---------------------------|------------------------|
| None | MSDS | None | None |

Secure Job Analysis

| Tasks | Unwanted Incidents | Precautions | Actions |
|---|---|--|---|
| Loading of bacterial cultures and chemicals | Spill of chemicals and bacterial culture from samples | Use specific tubes which tolerate proper chemicals and speed The use of gloves and lab coat are strongly recommended. | Use wipe paper to wipe off spill and ethanol 70% to decontaminate surfaces. And autoclaving of rotor in worst cases of contamination |
| Usage of Rotor | Expired age of rotor corrosion | The use of a LOG book or other means of determining rotor usage is required | LOG book |